




# LIQUID INSPECTION

Adelphi Manufacturing's Apollo Liquid Viewers offer a simple yet effective method of detecting particulate matter present in liquids.

## Apollo I Liquid Viewer

Designed for use in pharmaceutical and laboratory environments, the Apollo I Liquid Viewer has a polarised inspection light for the detection of particulate matter in transparent liquids contained in ampoules, vials and bottles. Particulate contamination can be unintentionally present in injections and intravenous infusions, and this unit provides a simple method of detection.

To operate, simply hold the container between the light panel and the polarising screen, shake or swirl the container, and any particles will show as bright specks against a dark background.



### Features:

- Compact and hygienic stainless steel body
- Adjustable viewing angle to suit individual working preferences
- Adjustable magnifier
- LED light source for optimum bulb life
- Tool free assembly
- Offers a clear, polarised, and magnified image
- Includes padded storage bag for safe storage & transport

| Item code         | Specification          |
|-------------------|------------------------|
| 6031300           | APOLLO I LIQUID VIEWER |
| Length (mm)       | 320                    |
| Width (mm)        | 295                    |
| Height (mm)       | 420 (max)              |
| Weight (kg)       | Nett: 4.9 Gross: 7.2   |
| Carton size (cm): | 50 x 41 x 27           |
| Electricity:      | 110/240V AC 50/60HZ    |




## Apollo II Liquid Viewer

Designed for use in pharmaceutical and laboratory environments, the Apollo II Liquid Viewer provides a simple method for the detection of visible particles in ampoules, vials, and bottles.

The unit is operated by simply swirling each individual container in front of both panels, and observing for any particulate contamination (other than air bubbles) unintentionally present in the solution. In front of the black panel, light coloured particulates will be visible, and in front of the white panel, dark coloured particulates will be visible.

This unit conforms to the European Pharmacopoeia Specification 2.9.20.



### Features:

- Adjustable lamp unit for comfortable working height
- Uses two energy efficient LED lamps
- Fold up frame for space saving & easy storage
- Removable diffuser for use with coloured glass
- Conforms to European Pharmacopoeia 2.9.20
- Intuitive assembly with no tools required

| Item code         | Specification           |
|-------------------|-------------------------|
| 6031510           | APOLLO II LIQUID VIEWER |
| Length (mm)       | 714                     |
| Width (mm)        | 420                     |
| Height (mm)       | 450                     |
| Weight (kg)       | Nett: 10 Gross: 7.9     |
| Carton size (cm): | 76 x 52.5 x 26.5        |
| Electricity:      | 100/240V AC 50/60HZ     |

We also offer an Apollo II High Intensity Specification unit, ideal for coloured glass containers, which uses an updated three LED bulbs, giving a light level of 8,000 – 10,000 lux.